---

---

**QuikRead®**

# Dispenser

**0,8 ml**

**REF** 06075

Cat.No.

**Operating Manual**

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**Mode d'emploi**

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**Istruzioni d'uso**

**Brugsanvisning**

**Bruksanvisning**

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
**Navodila za uporabo**

**Orion Diagnostica**

---

---

### Components

1. Handle
2. Outer sleeve
3. Plunger button
4. Valve head (S 28)
5. Intake tube with securing nut
6. Discharge tube securing nut
7. Discharge tube
8. Pressure compensation channel

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---

---

## 1. Intended Use and Operating Limitations

The dispenser should be used for QuikRead buffers only.

When used correctly, the dispensed liquid comes into contact with the following materials:

Borosilicate glass, FEP, ETFE, PFA, PTFE, Platinum-Iridium.

## 2. Safety Instructions

This Manual does not purport to address every safety issue which may arise during use. It is the user's responsibility to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Read this manual carefully and follow step by step.
2. Observe QuikRead reagent instructions supplied with each reagent kit. Ensure that the assay needs 0.8 ml buffer solution and that the QuikRead 0.8 ml Dispenser is used.
3. Follow the general instructions for hazard prevention and safety regulations, e. g., wear protective clothing, eye protection and gloves.
4. Use the dispenser only for dispensing buffers as part of the QuikRead system. To avoid contamination of the buffers, use a separate dispenser for each type of QuikRead assay.
5. Never use force on the dispenser.
6. While dispensing, the discharge tube must always point away from the user or other persons. Avoid splashes. Only dispense into suitable vessels.
7. Inspect the dispenser and its discharge tube regularly to detect any leakage.
8. Clean the dispenser regularly to ensure its smooth operation.
9. Use only original parts and accessories from the manufacturer. Do not carry out any technical alterations.
10. Clean the instrument before disassembling.
11. In case of difficulty (e.g., piston moves stiffly, leakage), immediately stop dispensing. Follow the Troubleshooting (chapter 8).

## 3. Instrument specifications

Accuracy	$\leq \pm 3 \%$ (0.8 ml)
Precision	$\leq 2 \%$ (0.8 ml)

Accuracy and precision (based on nominal volume) were determined with distilled water at the reference temperature of 20 °C.

---

---

## 4. Assembling

Before assembling, carefully inspect your QuikRead 0.8 ml Dispenser for any shipping damage. Check that the contents are complete (see page 2).

Any defects or shortages should be reported to the supplier.

The dispenser is designed to fit QuikRead buffer bottles and has been preset at the proper buffer dispensing volume of 0.8 ml.

1. Observe Safety Regulations (see chapter 2) and Operating Limitations (see chapter 1).
2. Screw the dispenser onto the buffer bottle. Do not screw too hard.
3. Affix to the dispenser a label with the name of the QuikRead assay you are going to use. To avoid contamination use a separate dispenser for each type of QuikRead assay.

**Warning!**

Avoid splashing of reagent!

4. Place a QuikRead cuvette below the discharge tube.
5. The discharge tube must point away from the user at all times.

**Note!**

Before using the instrument for the first time, ensure it is rinsed carefully or discard first few samples dispensed.

## 5. Priming

The dispenser should be primed before use and after each cleaning.

1. Place a suitable vessel below the discharge tube.

**Warning!**

Reagent may splash!

2. Gently push down and release the plunger button several times.
3. Check that there are no air bubbles in the discharge tube.

Every QuikRead buffer bottle contains excess buffer to allow sufficient priming.

---

---

## 6. Dispensing

1. Place a QuikRead cuvette below the discharge tube.
2. Slowly steadily depress plunger button, release gently. No air bubbles should be visible during dispensing.
3. Confirm the volume dispensed into the cuvette: the buffer surface should be between the two lines marked on the cuvette.

### **Note!**

If the piston moves stiffly or is difficult to move, stop dispensing immediately and clean the instrument (see chapter 7).

## 7. Cleaning and maintenance

Dispenser will only function safely if cleaned regularly, once a week is recommended. Be sure to clean the dispenser:

1. Immediately, if the plunger button moves stiffly
2. Before starting on a new reagent kit
3. Before longer downtimes
4. Before any maintenance and repair work

### **Warning!**

Dispenser, intake tube and discharge tube may be filled with reagent. To prevent injury from chemicals, always wear eye protection, protective clothing and protective gloves. Always point the discharge tube outlet away from the user or other people. Reagent may splash!

### 7.1. Emptying

1. Unscrew the dispenser and lift it high enough so that the intake tube is no longer immersed in liquid.
2. Cautiously tap the intake tube against the inside of the bottle so that any remaining reagent runs out.
3. Place a suitable vessel below the discharge tube.
4. Empty the discharge tube by repeated dispensing movements.

## 7.2. Cleaning procedure

1. Mount the empty dispenser (see chapter 7.1.) on a bottle filled with water.
2. Place a suitable vessel below the discharge tube.
3. Rinse the dispenser by repeated dispensing movements.
4. Completely empty the dispenser (see chapter 7.1.).
5. Assembling and priming as described in chapter 4 and 5.

## 8. Troubleshooting

Trouble	Possible causes	Action to be taken
<b>Aspirating and/or Dispensing not possible</b>	Sticking valves	Clean the dispenser (see chapter 7)
<b>Air is aspirated / air bubbles appear in the discharge tube</b>	Intake or discharge tube not firmly connected  Wrong or damaged intake or discharge tube	Push on intake or discharge tube more firmly. Screw on securing nut more firmly.  Use only original parts
<b>Plunger button moves stiffly</b>	Formation of crystals	Clean the dispenser (see chapter 7)
<b>Dispensed volume is too low Vacuum inside buffer bottle</b>	Valves are leaking  Pressure compensation channel is blocked.  Dispensing too fast	Clean the dispenser (see chapter 7)  Check the pressure compensation channel and clean it if necessary.  Observe dispensing speed.

If the measures do not help, send your dispenser to your local supplier for service. Please give full details of the fault and liquids dispensed. Clean the dispenser according to the normal procedure before despatching it back.

---

---

## 9. Warranty

The manufacturer's warranty for QuikRead 0.8 ml Dispenser covers defects in materials or workmanship for a period of six months from the date of purchase.

The manufacturer agrees to repair or replace at his sole discretion the dispenser if it becomes inoperative because of a failure of any internal part. The warranty does not cover damage caused by use not in accordance with instructions.

The manufacturer is under no obligation to modify or update the dispenser once it has been manufactured, unless a manufacturing defect is identified.

In the event of dispenser malfunction, contact your local supplier to obtain service during the six-month limited warranty period.

The QuikRead 0.8 ml Dispenser is manufactured for Orion Diagnostica, Finland by VITLAB GmbH, Germany.

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0,8 ml Dispensor QuikRead proizvaja podjetje VITLAB GmbH, Nemčija za podjetje Orion Diagnostica, Finska.
Orion Diagnostica Oy