







151053-7

- English
- Deutsch
- Français
- Eesti

QuikRead go[®] iFOBT

Assay Procedure • Testverfahren

Measure		RESULT
8	FOB	Hb in buffer: > 1000 ng/ml Hb in faeces: > 200 µg/g
Patient ID:	XXXXXXXX	Measurement time: 2013-06-11 12:19
Test:	iFOBT	Result info
	Choose Result info to view result information. Remove cuvette to perform a new measurement.	
Exit	Print	New measurement

1 Intended use

QuikRead go iFOBT is a test for the detection and quantification of haemoglobin originating from occult blood in human faeces in cases of excess bleeding from the lower gastrointestinal tract seen e.g. in ulcerative colitis, Crohn's disease, polyps, adenomas and colorectal cancer (CRC). The test can be used in routine physical examinations and in screening programmes for CRC. The measurement is carried out with the QuikRead go instrument.

For *in vitro* diagnostic use only.

2 Summary and explanation of the test

The faecal haemoglobin (Hb) concentration varies with gender and age, being higher for males than females and increasing by age for both genders.¹ The European Group on Tumour Markers (EGTM) recommends the use of quantitative, immunochemical faecal occult blood tests, iFOBT (faecal immunochemical test, FIT), which enable the clinician to set the appropriate cut-off value.² It should also be possible to express faecal Hb results as micrograms of Hb per gram of faeces, especially when comparing various methods.³

A high FOB value may indicate that the patient has large polyps or colorectal cancer (CRC). It usually takes 10 to 15 years for the first abnormal cells to develop into cancer. If these are detected and removed in time, the patient has near 100% reduction in risk of developing CRC.⁴ The lifetime risk of developing colorectal cancer (CRC) is about 1 in 20.⁵ In Europe, the annual incidence of CRC is predicted to increase from 432 000 new cases in 2008 to over 500 000 in 2020. Annual mortality is estimated to increase in the same proportion, reaching approx. 250 000 cases.⁶ Several studies have shown that CRC incidence and mortality can be significantly reduced by systematic screening programmes using immunochemical FOB tests. Both the Council of Europe and the American College of Physicians recommend regular CRC screening for people with no known risk factors to begin at the age of 50.^{7,8} Immunochemical FOB tests – iFOBT (FIT) – have brought significant improvement in analytical specificity.⁸

3 Principles of the procedure

QuikRead go iFOBT is an immunoturbidimetric test based on microparticles coated with porcine anti-human haemoglobin antibodies. The haemoglobin present in the sample reacts with the microparticles, and the resultant change in the turbidity of the solution is measured by the QuikRead go instrument.

The faecal sample is first prediluted in the QuikRead iFOB Sampling vial. The prediluted sample is added to the prefilled QuikRead go iFOBT cuvette. Test calibration data is contained in the barcode label of the cuvette and read automatically by the QuikRead go instrument before the test is started.

The instrument will indicate the presence of haemoglobin by displaying the detected concentration in terms of ng/ml, µg/g of faeces or as text FOB POSITIVE, depending on the chosen application in the instrument.

4 Reagents

Kit contents

Component name and origin	Symbol	QuikRead go® iFOBT Cat. no. 151051, 50 tests
FOB Reagent Caps ORIG PIG	REAG CPS	2 x 25
Buffer in prefilled cuvettes	BUF	2 x 25 x 1 ml
Instructions for use		

Reagents contain sodium azide, please see Section 5 "Warnings and precautions".

Storage and stability of the reagents

Kit component	Storage at 2...8°C	Storage at 18...25°C
Reagent caps in unopened tubes	Until the expiry date of the reagent tube	24 hours per day – 1 month 8 hours per day – 3 months
Reagent caps after first opening	Until the expiry date of the reagent tube	24 hours per day – 1 month 8 hours per day – 3 months
Prefilled cuvettes in an unopened foil pouch	Until the expiry date of the cuvette pouch	Until the expiry date of the cuvette pouch
Unopened prefilled cuvettes after the first opening of the foil pouch	6 months	3 months
Opened prefilled cuvette	2 hours	2 hours

Mark the date of opening the foil pouch on the cuvette rack.

Reagent preparation and storage conditions

All reagents are ready to use. Keep the FOB Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

Reagent deterioration

The product should be used only if the buffer volume in the cuvette is correct. Check that the liquid surface is between the two lines marked on the cuvette. Do not use a cuvette with visible dirt in the buffer.

5 Warnings and precautions

Health and safety information

- For *in vitro* diagnostic use only.
- Do not smoke, eat or drink in premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of water.
- Do not drink the liquid in the sampling tube.
- All patient samples and controls should be handled as potentially infectious material.
- Lyophilised reagents contain < 1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501).
Reconstituted and liquid reagents contain < 0.1% sodium azide, which is not considered a harmful concentration.
Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Disposal: see Chapter 14.

Analytical precautions

- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.

- The QuikRead go iFOBT kit is intended solely for use with the QuikRead go instrument.
- Do not mix reagents from different lot numbers or different tests. The components are disposable; never reuse components already used for performing a test.
- When opening a kit for the first time, ensure that the foil pouches (2 pcs) protecting the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes inside it. In addition, always ensure before using an individual cuvette that its foil covering is intact.
- Do not touch the clear flat surfaces on the lower part of the cuvette (optical part). Discard any cuvettes with fingerprints.
- QuikRead FOB Reagent Caps are colour-coded brown to distinguish them from other QuikRead analytes.
- Do not splash liquid into the measuring well of the instrument.

6 Sample collection and preparation

The patient does not have to undertake any special diet prior to faeces sampling. Provide the patient with the iFOB sampling tube, a patient label and a copy of the following sampling instructions. Explain the instructions to the patient.

Do not drink the liquid in the sampling tube. Take care not to bend or break the tip of the tube.

- Obtain a sample of faeces in a clean container or on clean collecting paper. Avoid contact with the toilet water when you are collecting the sample. If there is visible blood in the faeces sample or in case of diarrhea, do not collect a sample but wait for the next bowel movement and collect a sample from that.
- Unscrew the top of the sampling tube containing the sample collecting stick. Collect the sample by poking the sample collecting stick (attached to the tube cap) in three (3) different places on the faeces sample.
- Insert the sample collecting stick back into the sampling tube and screw the cap on firmly.
- Write your name and the date on a separate label and attach the label to the sampling tube.
- Return the iFOB sampling tube to the laboratory or doctor as soon as possible after sampling. Protect the sampling tube from breakage or leakage during transportation.
- The sample in the sampling tube can be stored in a refrigerator (2...8°C) for up to 7 days or at room temperature (max. 25°C) for up to 5 days. Storage in refrigerator is however preferred.

Allow the sample to reach room temperature (18...25°C) before testing.

7 Procedure

Materials required but not supplied

	Cat. no.
QuikRead go® Instrument	133893
QuikRead® iFOB Sampling Set (50 sampling tubes)	151052
QuikRead® FOB Positive Control	06027
QuikRead go® iFOBT Control Quantitative	152390

Assay procedure

Open the foil pouch protecting the cuvette rack and mark the opening date on the cuvette rack.

Sample handling (Fig. 1–5)

Process the sample immediately after receiving it. The sampling tube contains formaldehyde (may cause cancer and may produce an allergic reaction). For full instructions, see QuikRead iFOB Sampling set instructions for use.

- 1 A prefilled cuvette must reach room temperature (18...25°C) before use. This will take 15 minutes for an individual refrigerated (2...8°C) prefilled cuvette. Do not touch the clear flat surfaces in the lower part of the cuvette (the optical part). Remove the foil cover of the cuvette. Beware of splashing the liquid. Condensed water on the foil surface has no effect on results. The test must be performed within two (2) hours from opening the cuvette. The liquid surface should lie between the two lines marked on the cuvette.
- 2 Check that each sampling tube is marked with the patient's ID and sampling date. If the label is missing, mark the tube with the patient's ID and sampling date. Shake the sampling tube and break the tip of it. Foam may form during mixing. Beware of splashes when the tip is broken. Wipe the foam off the tip.
- 3 Discard the first three (3) drops.
- 4 Then dispense three (3) drops into the cuvette.
- 5 Close the cuvette tightly with the FOB Reagent Cap. Do not press down the brown inner part of the reagent cap. After the sample has been added to the buffer, the solution is stable for two (2) hours. Keep the cuvette in an upright position, do not shake.

Analysing the sample (Fig. 6–8)

For full instructions, see QuikRead go instrument instructions for use. The display will guide you through the test.

- 6 On the QuikRead go instrument display, select **Measure**.
- 7 Place the cuvette in the instrument with the barcode facing you, as shown in figure 7. The instrument recognises the sample type. The display shows how the measuring progresses. The instrument first measures the sample blank and the Hb concentration within two (2) minutes.
- 8 The result appears on the display when the measurement is completed and the cuvette automatically rises from the measurement well.

Analysing the control

Regular use of FOB Positive Control (Cat. no. 06027) and/or QuikRead go iFOBT Control quantitative (Cat. No. 152390) is recommended.

Precise instructions on the use of controls are provided in the package insert of the control vials. Analyse the control samples as patient samples but select **Quality Control** on the instrument display. The result will be saved as control measurement.

The result displayed depends on the instrument settings.

Concentration of haemoglobin in the sampling tube buffer	Quantitative result			Qualitative result
	Depending on chosen instrument settings, the display shows both values or only either of the values			
< 50 ng Hb/ml buffer	< 50 ng Hb/ml buffer	and/or	< 10 µg Hb/g faeces	< 15 µg Hb/g FOB NEGATIVE
50–1000 ng Hb/ml buffer	50–1000 ng Hb/ml buffer	and/or	10–200 µg Hb/g faeces	
> 1000 ng Hb/ml buffer	> 1000 ng Hb/ml buffer	and/or	> 200 µg Hb/g faeces	≥ 15 µg Hb/g FOB POSITIVE

8 Quality control

Regular use of FOB Positive Control (Cat. no. 06027) and/or QuikRead go iFOBT Control quantitative (Cat. no. 152390) is recommended.

QuikRead FOB Positive Control (Cat. no. 06027):

- When QuikRead go instrument result settings are chosen to be qualitative, "FOB POSITIVE" should be shown on the display.
- When QuikRead go instrument result settings are chosen to be quantitative, "Hb in buffer > 1000 ng/ml and/or > 200 µg Hb/g faeces" should be shown on the display.

QuikRead go iFOBT Control quantitative (Cat. no. 152390):

- When QuikRead go instrument result settings are chosen to be qualitative, "FOB POSITIVE" should be shown on the display.
- When QuikRead go instrument result settings are chosen to be quantitative, a result within the range indicated on the control vial label should be displayed.

9 Results' interpretation

The instrument measures the concentration of Hb in the cuvette buffer and calculates the amount of Hb per gram of faeces based on the assumption that the sampling has been done correctly according to the instructions for use (~10 mg faeces in 2 ml buffer). If the amount of faeces in the sampling tube is higher or lower than instructed, the Hb value per g faeces will deviate in the same proportion.

10 Limitations of the procedure

QuikRead go iFOBT measures haemoglobin in the faecal sample. Since bowel lesions may not bleed at all or may bleed intermittently, or the blood may be unevenly distributed in the faecal sample, a test result may be negative even in cases with occasional bleeding. Healthy persons may have blood in their faeces because of physical stress or intake of certain medications that may cause gastrointestinal irritation resulting in gastrointestinal bleeding. As with any faecal occult blood test, QuikRead go iFOBT may not be considered conclusive proof of gastrointestinal bleeding or pathology. The faecal occult blood test results should always be interpreted in the light of the patient's symptoms and current diagnostic criteria. The faecal occult blood tests are not intended to replace other diagnostic procedures such as sigmoidoscopy, colonoscopy or X-ray studies but are to be regarded as preliminary screening methods or aids to diagnosis.

11 Expected values

In a healthy person, there should be no haemoglobin or a minimal amount of haemoglobin (physiological bleeding) in faeces, see "Limitations of method". The result obtained should be <75 ng Hb/ml buffer, <15 µg Hb/g faeces or FOB NEGATIVE.

12 Performance characteristics**Precision and reproducibility**

A precision study was performed according to Clinical and Laboratory Standards Institute (CLSI) guideline EP05-A3:2014 using the single-site precision evaluation protocol. Samples containing three different levels of human hemoglobin in feces suspension were evaluated. The test was conducted during 20 testing days, two runs per day, and two replicate measurements per run for each sample level using one reagent lot, one instrument and one user. Samples were dispensed using three drops from the sampling tube and a control was dispensed using a pipette.

Sample	Mean Hb (ng/ml)	N	Repeatability		Between run		Between day		Within laboratory	
			SD	CV%	SD	CV%	SD	CV%	SD	CV%
Sample 1	64	80	1.7	3.0	5.9	9.2	0.0	0.0	6.2	9.7
Sample 2	269	80	4.2	1.5	17.5	6.5	8.7	3.2	20.0	7.4
Sample 3	712	80	65.5	9.2	88.9	12.5	0.0	0.0	110.4	15.5
Control	177	80	3.0	1.7	1.1	0.6	0.0	0.0	3.2	1.8

Abbreviations: CV%, coefficient of variation expressed as a percentage; SD, standard deviation; N, number of results.

Linearity

Linearity was determined according to Clinical and Laboratory Standards Institute (CLSI) guideline EP06-A:2003. For quantification of Hb by QuikRead go iFOBT 151051 the method has been demonstrated to be linear across the measuring range within allowable nonlinearity of ±20 %.

Limit of Detection (LoD)

The LoD is 2,5 µg Hb/g faeces, determined according to CLSI guideline EP17-A2:2012, based on 300 determinations, with 150 blank and 150 low level replicates; and a LoB of 1,2 µg Hb/g faeces.

Limit of Quantitation (LoQ)

The LoQ is 9,5 µg Hb/g faeces, determined according to CLSI guideline EP17-A2:2012 based on 90 determinations and a TE goal of 20% using the Westgard model.

Method comparison

A method comparison was performed between another commercially available immunological FOB method and the QuikRead go iFOBT method. Human patient faeces samples were analysed. A summary of obtained qualitative results is presented in the table below.⁹

Test method		Comparison method	
		+	-
QuikRead go iFOBT	+	40	6
	-	1	65
Sensitivity to comparison method		97.6%	
Specificity to comparison method		91.5%	
Agreement		93.8%	

Antigen excess

Haemoglobin concentrations less than 500 µg/ml (<100 mg haemoglobin in 1 g faeces) do not give falsely negative results.⁹

Interfering substances

No interference with vitamin C at concentrations less than 100 µg/ml, an amount corresponding to a daily intake of approx. 2 g ascorbic acid, has been observed. Iron supplementation does not interfere with the measurement.⁹

Measuring range

The QuikRead go instrument with software version 6.1.5 or higher will give quantitative or qualitative results depending on the QuikRead go instrument's settings.

For the quantitative test, the Hb measuring range is 50–1000 ng/ml test solution and/or 10–200 µg/g faeces. In cases when the result is above the measuring range, the result will be displayed as > 1000 ng/ml and/or > 200 µg/g faeces. If the result is below the measuring range, the result will be displayed as < 50 ng/ml and/or < 10 µg/g.

For the qualitative result the cut-off value is 75 ng/ml test solution, which corresponds to 15 µg/g faeces. In cases when the result is at the cut-off value or above the cut-off value, the result will be displayed as "FOB POSITIVE" and in cases when the result is below the cut-off value, the result will be displayed as "FOB NEGATIVE".

13 Traceability

Each batch is tested against reference material of human origin. The haemoglobin content of this reference material has been verified with a method comparable to the International Council for Standardization in Haematology (ICSH) reference method.

14 Disposal

- Dispose of contents according to national and local law.
- All patient samples, sampling devices, controls, used tubes, caps and cuvettes should be handled and disposed of as potentially infectious material.
- Materials of the components of the QuikRead 101 and QuikRead go kits:
 - Paper: Instructions for use
 - Cardboard: Kit box including its inner parts
 - Plastic: Cuvettes, reagent caps, foil covering cuvette rack, cuvette rack, plungers, extraction vials and tubes, plunger and capillary tubes, swabs
 - Glass: Capillaries
 - Metal: Reagent cap tubes, cuvette lids, plunger and capillary tube caps
- Not to be recycled: Lids of reagent cap tubes (several) and magnetic cards (PVC)
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

15 Troubleshooting

Error messages

Error messages displayed by the QuikRead go instrument are listed below. For further details of error messages, consult the QuikRead go instrument instructions for use.

Error message	Corrective action
Cuvette position not correct. Remove cuvette.	Remove the cuvette and re-insert it into right position.
Measurement prohibited. Please check the reagent cap.	Check that the cuvette has a reagent cap and the inner brown part of the cap is not pressed down.
Measurement prohibited. Kit lot expired.	Discard expired kit lot. Use a new one.
Measurement prohibited. Cuvette temperature too low.	Let the cuvette warm up to room temperature (18...25°C). Test the same cuvette again.
Measurement prohibited. Cuvette temperature too high.	Let the cuvette cool down to room temperature (18...25°C). Test the same cuvette again.
Test cancelled. Blank too high.	Test the same cuvette again. The blanking process has not been completed or the sample may contain interfering substances. In the latter case, the test cannot be completed.
Test cancelled. Unstable blank.	
Test cancelled. Error in reagent addition.	Perform a new test. There has been some problem during reagent addition. Make sure that the cap is properly closed.
Test cancelled. Instrument failure.	Perform a new test. If this message appears often, contact your local distributor or call customer service.
Error code xx. Please restart the instrument.	Restart the instrument. If the error message appears again, contact your local distributor or call customer service and give the error code.
Error code xx. Please call customer service.	Contact your local distributor or call customer service and give the error code.

Unexpected low and high results

Possible reasons for unexpected low and high results are listed in the table below.

Problem	Possible cause	Corrective action
Unexpected low result	The sample volume in the cuvette is too small.	Perform a new test. Make sure that the correct number of drops (3) is added.
	The amount of sample taken in sampling is too small.	Perform a new test. Make sure that the sampling instructions are followed (see section 6: Sample collection and preparation).
Unexpected high result	The sample volume in the cuvette is too large	Perform a new test. Make sure that the correct number of drops (3) is added.
	The sample has been added to the sampling tube more than once.	Perform a new test. Make sure that the sampling instructions are followed (see section 6: Sample collection and preparation).
	The buffer volume is too low.	Perform a new test. Make sure that the volume is correct by checking that the liquid surface is between the two lines marked on the cuvette.
	The cuvette is dirty.	Perform a new test. Do not touch the clear flat surfaces in the lower part of the cuvette.
Unexpected low/high result	Components from different kit lots or tests were used.	Perform a new test. Make sure that all the reagents are from the same kit lot.
	Incorrect reagent storage.	Perform a new test. Make sure that the reagents are stored according to the instructions for use.

References • Referenzen

1. McDonald PJ et al. Faecal hemoglobin concentrations by gender and age: implications for population-based screening for colorectal cancer. Clin Chem Lab Med. 2012; 50 (5): 935-40
2. Duffy MJ et al. Use of faecal markers in screening for colorectal neoplasia: a European group on tumor markers position paper. Int J Cancer. 2011; 128 (1): 3-11
3. Fraser CG et al. A Proposal to Standardize reporting Units for Fecal Immunochemical tests for Hemoglobin. J Natl Canc Inst. 2012; 104 (11): 810-4
4. Quirke P et al. Quality assurance in pathology in colorectal cancer screening and diagnosis – European recommendations. Virchows Arch. 2011; 458 (1): 1-19
5. van Roon AH et al. Are Fecal Immunochemical Test Characteristics Influenced by Sample Return Time? A Population-Based Colorectal Cancer Screening Trial. Am J Gastroenterol. 2012; 107 (1): 99-107
6. Qaseem A et al. Screening for colorectal cancer: a guidance statement from the American College of Physicians. Ann Inter Med. 2012; 156 (5): 378-86
7. Globocan (<http://globocan.iarc.fr/>)
8. European Commission: European Guidelines for Quality Assurance in Colorectal Cancer Screening and Diagnosis – First Edition. Segnan N, Patrick J, von Karsa L (eds), 2010. Luxembourg: Publications Office of the European Union
9. Data on file, Aidian Oy

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QuikRead go[®] is a registered trademark of Aidian Oy.

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