## QuikRead go<sup>®</sup> wrCRP Control, Cat. No. 153038

153069-1

*English*

<p><b>Intended use</b></p> <p>Orion Diagnostica QuikRead go<sup>®</sup> wrCRP Control is intended for quality control of wrCRP assays by the QuikRead go<sup>®</sup> Instrument.</p> <p><b>For <i>in vitro</i> diagnostic use.</b></p> <p><b>Principles of the procedure</b></p> <p>Use the control in the same way as unknown <b>plasma/serum</b> specimen would be used, following the instructions for the assay procedure.</p> <p><b>Reagents</b></p> <p>The control is human-based material. It should be handled in the same way as patient sample.</p> <p><b>Preparation</b></p> <p>QuikRead go wrCRP Control is ready-to-use. Allow the material to stabilise for 30 minutes at room temperature (18...25°C). Mix well. Do not shake.</p> <p><b>Storage</b></p> <p>Stable until the expiry date when stored at 2...8°C. Opened vial: Stable for 2 months at 2...8°C. Do not use any vial that has visible evidence of microbial growth.</p> <p><b>Expected values</b></p> <p>The individual result should be within the range indicated on the control vial label. Each laboratory should establish its own mean and precision data.</p>	
<p><b>Warnings and precautions</b></p> <ul style="list-style-type: none"><li>Reagents contain &lt;0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.</li> <li>Avoid using the last drops of control because opening of the vial many times may result in evaporation, which may bias the results. If the control is not in the acceptable range, take a new bottle.</li> <li>Do not use the product after the expiry date marked on the outer package.</li> <li>Do not exceed the stability periods for opened reagents.</li> <li>QuikRead go wrCRP Control contains material of human origin, which has been tested and found negative for human immunodeficiency virus antibodies (HIV 1 and 2), hepatitis C virus (HCV) and hepatitis B surface antigen (HBsAg). As a precaution, handle QuikRead go wrCRP Control as if capable of transmitting infectious agents.</li></ul>	
<p><b>Disposal</b></p> <ul style="list-style-type: none"><li>Dispose of contents according to national and local law.</li> <li>Materials of the components</li> <li>Paper: Instructions for use</li> <li>Cardboard: Outer box</li> <li>Glass: Control vials</li> <li>Plastic: Vial caps</li> <li>Rubber: Stopper</li> <li>When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.</li></ul>	

<p><b>QuikRead go<sup>®</sup> wrCRP Control, Cat. No. 153038</b></p> <p>QuikRead go<sup>®</sup> wrCRP Kontrolle, Kat. -Nr. 153038</p>	
<p><i>Deutsch</i></p>	

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**Warnings and precautions**

- Reagents contain <0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Avoid using the last drops of control because opening of the vial many times may result in evaporation, which may bias the results. If the control is not in the acceptable range, take a new bottle.
- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- QuikRead go wrCRP Control contains material of human origin, which has been tested and found negative for human immunodeficiency virus antibodies (HIV 1 and 2), hepatitis C virus (HCV) and hepatitis B surface antigen (HBsAg). As a precaution, handle QuikRead go wrCRP Control as if capable of transmitting infectious agents.

**Disposal**

- Dispose of contents according to national and local law.
- Materials of the components
- Paper: Instructions for use
- Cardboard: Outer box
- Glass: Control vials
- Plastic: Vial caps
- Rubber: Stopper
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

<p><b>QuikRead go<sup>®</sup> wrCRP Control, Cat. No. 153038</b></p> <p>QuikRead go<sup>®</sup> wrCRP Kontrolle, Kat. -Nr. 153038</p>	
<p><i>Deutsch</i></p>	

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