## Intended use

Hygicult CF slides are intended for presumptive detection of coliform bacteria. The slide is covered on both sides with VRB (Violet Red Bile) Agar which allows coliforms to grow as bright red colonies.

The test can be performed on-site for monitoring different types of materials, both solid and liquid. As required, the slides can be used as convenient transport media for samples.

**Note:** The limit values for microbial count in normal drinking water are too low to be detected by the Hygicult method.

## Contents of the kit

<b>Hygicult CF</b>	<b>Cat. No. 06036</b>
Test slides	10 pcs
Labels	10 pcs
Instructions for use	1 pc

## Typical formulation

<b>VRB Agar</b>	
Peptone	Neutral red
Yeast extract	Crystal violet
Sodium chloride	Agar agar
Lactose	Water
Bile salts	

## Warnings and precautions

Do not use product beyond the expiry date marked on the kit. Do not use the kit if you notice

- discoloration or dehydration of the growth medium
- detachment of the growth media from the plastic slide
- evidence of bacterial or fungal growth

Do not touch the growth because any colony growing on the slide may be pathogenic.

## Storage

Store the kit at room temperature (18...25°C) protected from draught, temperature fluctuations and light sources. Avoid storage near heat-generating appliances. Do not allow to freeze. The expiry date (year-month-date) is marked on the box and on the cap of each slide.

## Sampling

To avoid contamination, the growth medium should not come into contact with any other material than the one to be tested. On the other hand, it is important that the growth medium makes full contact with the material to be tested. After sampling screw the slide tightly back into the tube.

## Contact inoculation (Fig. 1a, 1b)

Solid surfaces can be tested by pressing each side of the slide firmly against the surface for three or four seconds. The slide should be held still during pressing. The hinged design offers ease of use.

## Dipping (Fig. 2)

Fluid samples are tested by dipping the slide in the liquid for three or four seconds. Blot the last drops on absorbent paper.

## Swabbing (Fig. 3)

Semisolid materials or objects that are difficult to reach can be tested by carefully rolling a sterile swab over an area delimited using e.g. a frame. If the object is dry, the swab should first be moistened with sterile water. The moistened swab can also be used for obtaining samples from powders (e.g. spices) or viscous fluids.

After swabbing the sample area, roll the swab gently over the agar surfaces of the slide from left to right and from bottom to top.

## Incubation (Fig. 4)

Incubate the slide tightly enclosed in its tube at 35...37°C for 24–48 hours.

## Interpretation of results (Fig. 5)

Remove the slide from its tube after incubation and determine the microbial count (number of colony forming units, CFU) and examine the colour reactions by comparing with the model chart.

Coliform bacteria grow on VRB Agar as bright red colonies. The growth of gram-positive organisms is inhibited.

The following levels can be considered as a rough basis for evaluating the degree of contamination.

	<b>Contact inoculation</b>
Clean	0 CFU/slide
Contaminated	1–10 CFU/slide
Very contaminated	> 10 CFU/slide

The presence of coliform bacteria in cooked food always indicates mishandling of the product or inadequate hygiene.

## Limitations of the method

When used as a contact slide, Hygicult CF equals the contact plate method in sensitivity, whereas the dip and swab procedures have a detection limit of 1000 CFU/ml. The allowed total microbial concentration of normal drinking water is too low to be reliably detected using Hygicult CF.

Results obtained with different inoculation systems should not be compared. Valid comparisons can only be made among results obtained using the same technique on the same type of material.

## Disposal

Any growth on slides may be pathogenic. Used slides must therefore be disposed of by burning, autoclaving or immersion in a disinfectant overnight, always following local laws and regulations.

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Hygicult® is a registered trademark of Orion Diagnostica Oy.