Intended use

A dip slide culture method for diagnosing urinary tract infections (UTI). Uricult Vet is intended for veterinary use only.

Principles of the procedure

The Uricult Vet dipslide system is based on two culture media for detection of microbes causing urinary tract infections in animals. One side of the plastic slide is covered with green CLED medium and the other side with reddish-brown MacConkey medium. The CLED medium is intended for determining the total microbial count. The MacConkey medium is intended for detecting gram-negative microbes. This medium contains bile salts that prevent the growth of gram-positive microbes other than enterococci, which may grow as pinpoint colonies.

Reagents

Contents

| Uricult Vet | Cat. No. 06064 |
|-----------------------|----------------|
| Uricult Vet dipslides | 10 |
| Patient labels | 10 |
| Instructions for use | 1 |

Typical formulation

| CLED medium | | MacConkey medium | |
|-----------------|----------|------------------|-----------|
| Peptone | 10.0 g/l | Peptone | 20.0 g/l |
| Meat extract | 3.0 g/l | Lactose | 10.0 g/l |
| Lactose | 10.0 g/l | Neutral red | 0.075 g/l |
| L-Cystine | 0.13 g/l | Bile salts | 0.8 g/l |
| Bromthymol blue | 0.03 g/l | | |

Storage

Store Uricult Vet at 7...25°C protected from draught, temperature fluctuations and light sources. Avoid storage near heat-generating appliances. **Do not allow to freeze.** The expiry date is marked on the box.

Warnings and precautions

Uricult Vet is for **veterinary diagnostics** only.

Do not use the product beyond the expiry date marked on the box. Wear protective clothing and disposable gloves while handling samples or tests, and wash hands thoroughly afterwards. Do not use the Uricult Vet if you detect discoloration or dehydration of the agar, separation of the growth media from the plastic slide or evidence of bacterial or fungal growth. Because any colonies growing on Uricult Vet are actual or potential pathogens, do not touch the growth.

Sample collection and preparation

Ideally, urine for bacterial culture should remain in the bladder for four hours prior to sampling. The veterinarian will provide detailed instructions for "home sampling". The veterinarian can take a sample via catheterisation or cystocentesis, as required. Cystocentesis is preferred method when taking sample from animals suspected of having a UTI^{1,2}.

The urine sample may also be taken by holding a Uricult Vet slide directly in the animal's urine stream. **Please note that this method is not preferred.** In this method, it is important to ensure that the surfaces of both culture media become completely wet and that no contamination occurs from animal hair or other objects in conjunction with sampling.

The sample should be inoculated onto the Uricult Vet dipslide immediately after collection. The dipslide should then at once be returned into its protective tube and the cap closed.

If the urine sample needs to be stored prior to inoculation onto Uricult Vet, it should be kept refrigerated at 2...8°C no longer than 24 hours.

Anti-inflammatory medications may affect the result of the Uricult Vet test. Therefore, the test should not be performed until 48 hours after the final dose of medication.

Inoculated dipslide may be

- incubated immediately or
- stored at 2...8°C for up to 48 hours or
- transported to a laboratory for incubation and/or interpretation

Procedure

- Unscrew the slide from the tube without touching the agar surfaces.
- Holding Uricult Vet by the cap, dip the slide into freshly voided midstream urine so that the agar surfaces are totally immersed (a). If the volume of urine is too small for this or the sample has been taken by cystocentesis, the agar surfaces can be wetted by pouring urine (b) or placing drops from the syringe on them, followed by tilting to ensure complete wetting (c). Uricult Vet may also be inoculated directly in the animal's urine stream. Do this by removing a Uricult Vet slide from its protective tube and holding it by the cap in the urine stream. Tilt the slide from side to side to ensure that the surfaces of both culture media become wet while taking care that animal hair or other objects do not touch the surfaces of the culture media.
- Allow excess urine to drain from the slide (a) and blot the last drops on absorbent paper (b).
- Screw the slide tightly back into the tube.
- Fill in the patient label and attach it to the tube.
- Place the tube upright in an incubator (36±2°C) for 16–24 hours. The tube may also be sent to a laboratory for incubation.
- To obtain a colony forming units per millilitre (CFU/ml), remove the slide from the tube and compare the colony density with the model chart provided in the kit.

Note:

- Negative cultures and complicated UTI samples are recommended to always be incubated for an additional 24 hours to ensure that slow-growing bacteria are detected.
- The inoculated slide may be incubated immediately or stored or transported to a laboratory for incubation and interpretation. Storage or transportation should not exceed 48 hours at 7...25°C, after which Uricult Vet should be incubated at 36±2°C for 16–24 hours. If the slide has been

stored or transported for up to 48 hours, only the presence of growth and the colony count should be recorded from it; the colour reaction may be atypical.

- The inoculated slide may be incubated at room temperature for 24–72 hours, after which positive cultures may be interpreted.

Quality control

Quality control tests are performed on each lot of Uricult Vet dipslides at the time of manufacture. Good laboratory practice recommends the use of quality control. The following procedure is recommended:

- Prepare a 10⁵–10⁹ bacteria/ml suspension of each of the following bacteria in sterile saline:
 - Staphylococcus aureus* ATCC 25923
 - Escherichia coli* ATCC 25922
 - Proteus mirabilis* ATCC 12453
- Use the suspensions to inoculate the Uricult Vet dipslides, using the normal procedure.
- Interpret the results after 16–48 hour incubation as follows:

S. aureus ATCC 25923: Growth of colonies on the CLED medium only. Colonies ferment lactose, as indicated by the yellow colour of the colonies and the shift towards yellow of the medium.

E. coli ATCC 25922: Growth of yellow colonies with a shift towards yellow of the CLED medium and growth of pink-red colonies on the MacConkey medium.

P. mirabilis ATCC 12453: Growth of translucent colonies with a shift towards blue of the CLED medium and growth of colourless colonies on the MacConkey medium.

Results' interpretation

After incubation of the inoculated dipslide, the presence of bacteria is evidenced by colonies on the culture medium. Because a colony is the result of the multiplication of a single bacterial cell, the number of colonies indicates the number of colony forming units (CFUs) in the urine sample. The colony count is determined from the originally green CLED medium by matching the colony density with the model chart its most closely resembles. It is important to compare the number of colonies, not their size.

The low electrolyte concentration of the CLED medium prevents spreading of *Proteus* strains. Bromthymol blue and lactose in the medium allow the detection of lactose fermenting bacteria. Such lactose fermenting strains grow as yellow colonies and turn the culture medium yellow, whereas non-lactose fermenting strains are translucent with no colour change of the medium.

The originally brownish-red, selective MacConkey medium supports the growth of gram-negative bacteria, but certain enterococci may grow as pinpoint colonies on the medium³. The selectivity is accomplished by bile salts. Lactose fermenting bacteria grow as red and non-lactose fermenting bacteria as translucent colonies on the medium.

When the urinary bacterial content is high (≥ 10⁷ CFU/ml), the agar surfaces may become totally covered by confluent growth. This can be misinterpreted as a negative result. Therefore, any surfaces that appear negative should be examined under a reflecting light; absence of reflection indicates confluent growth. A bright light also allows very small colonies to be detected.

A mixture of different bacterial strains on the Uricult Vet is most likely due to contamination of the urine sample. Samples obtained by cystocentesis most often yield pure cultures whereas the voided samples are often contaminated, yielding mixed cultures⁴.

Limitations of procedure

Uricult is capable of detecting bacterial concentrations between 10³ and 10⁷ CFU/ml. The model chart allows the determination of colony counts to the nearest power of 10. When the chart is used according to instructions, colony counts show a 99% correlation with the conventional pour plate method⁵.

Expected values

High bacterial number in a properly collected and cultured sample indicates bacterial UTI¹.

| Method of sampling | Significant colony count CFU/ml ⁶ | |
|--------------------|--|-------------------|
| | Dog | Cat |
| Cystocentesis | ≥ 10 ³ | ≥ 10 ³ |
| Catheterisation | ≥ 10 ⁴ | ≥ 10 ³ |
| Voided urine | ≥ 10 ⁵ | ≥ 10 ⁴ |

Note that contamination of voided samples may cause growth of ≥ 10⁴ CFU/ml¹. Therefore voided samples should not be routinely used.

Performance characteristics


Uricult Vet • CLED medium⁷

| | | |
|-------------------|------|--|
| Number of samples | 140 | Reference method: Pour plate (Nutrient agar) |
| Sensitivity | 100% | |
| Specificity | 99% | |
| PPV | 98% | |
| NPV | 100% | |


Disposal

- Dispose of contents according to national and local law.
- All patient samples and used components should be handled and disposed of as potentially infectious material.
- Materials of the components:
 - Paper: Instructions for use, patient labels
 - Cardboard: Kit box
 - Plastic: Tubes, caps and dipslides
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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Model Density Chart • Mallitaulu • Modelldensitetsdiagram

The colony count is determined from the originally green CLED medium by matching the colony density with the model chart its most closely resembles. **It is important to compare the number of colonies, not their size.** A growth consisting of several species of bacteria is termed mixed flora and is most likely due to contamination of the urine sample.

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References • Viitteet • Referenser

1. Bartges JW. Diagnosis of urinary tract infections. *Vet Clin Small Anim.* 2004; 34: 823-933.
2. Weese JS, Blondeau JM, Boothe D *et al.* Antimicrobial use guidelines for treatment of urinary tract disease in dogs and cats: antimicrobial guidelines working group of the international society for companion animal infectious diseases. *Veterinary Medicine International* 2011; article ID: 263768. doi:10.4061/2011/263768
3. Granato PA. Evaluation of a dip slide device for enumeration of bacteria in urine. *Laboratory Medicine* 1980; 11: 4: 246-250.
4. Duijkeren E van, Laar P van, & Houwers DJ. Cystocentesis is essential for reliable diagnosis of urinary tract infections in cats *Tijdschr Diergeneeskd* 2004; 129: 394-396.
5. McAllister TA, Arnell GC, Barr W, & Kay P. Assessment of plain dipslide quantitation of bacteriuria. *Nephron* 1973; 11: 111-122.
6. Lulich JP, Osborne CA. Bacterial urinary tract infection. In: Ettinger SJ, Feldman EC, editors. *Textbook of veterinary internal medicine*. 4th edition. Philadelphia: WB Saunders, 1999: 1775-1788.
7. Arnell GC, McAllister TA, & Kay P. Detection of bacteriuria at room temperature. *Lancet* 1970; January 17: 119-121.

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Uricult® is a registered trademark of Orion Diagnostica Oy.

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