

QuikRead go[®] Connectivity

Take full advantage of connecting your QuikRead go[®] system

Take full advantage of connecting your QuikRead go[®] system

Saving patient results in an electronic format increases patient safety by removing risks of false manual recordings. QuikRead go test results transferred to a data handling system coupled with patient and/or operator identification allow to achieve this goal. QuikRead go can be connected to Laboratory Information Systems, Hospital Information Systems or to Middlewares using a bidirectional POCT1-A2 connection.

Because...

ENHANCED PATIENT SAFETY

- » Send results automatically to reduce manual errors
- » Download patient lists to reduce errors
- » Add comments with results

REMOTE INSTRUMENT MANAGEMENT


- » Lock/unlock QuikRead go instruments
- » Adjust instrument settings
- » Set QC limits
- » Update instrument software
- » Use quality control checking

REMOTE USER MANAGEMENT

- » Manage operator lists
- » Send information to QuikRead go operators

ALSO AVAILABLE

- » Unidirectional connection by LIS01-A2
- » Wireless connection with a WLAN USB stick



QuikRead go[®]



More information: www.oriondiagnostica.com
or contact Orion Diagnostica.

Product	Cat. no.
QuikRead go Instrument	133893
Printer for QuikRead go instrument	136203
Barcode reader for QuikRead go instrument	136204
QuikRead go WLAN adapter	148045



Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission.

QuikRead go[®] is a registered trademark of Orion Diagnostica Oy.

Orion Diagnostica Oy
Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland
Tel. +358 10 4261, Fax +358 10 426 2794
www.oriondiagnostica.com, www.quikread.com

8027-02GB, 01/2016. International version. Not intended for US market.