



Uricult® Trio

Model Chart



Uricult® Products

For diagnosis of urinary tract infections by demonstrating microbes in urine

Uricult 67404

- » CLED medium for total bacterial counts
- » Selective MacConkey medium for gram-negative bacteria



Uricult Plus 67465

- » CLED medium for total bacterial counts
- » Selective MacConkey medium for gram-negative bacteria
- » Selective enterococcus medium for enterococci, group B streptococci may also grow on this medium



Uricult Trio 68197

- » CLED medium for total bacterial counts
- » Selective MacConkey medium for gram-negative bacteria
- » E. coli medium for gram-negative, β -glucuronidase -producing organisms, such as *E. coli*. Certain yeasts may also grow on this medium.



| Product list | Cat. No. |
|----------------------------------|----------------|
| Uricult 10 tests | 67404 |
| Uricult Plus 10 tests | 67465 |
| Uricult Trio 10 tests | 68197 |
| Cultura M 230 V | 129997 / 05866 |
| Cultura M 110/115 V | 132754 / 06266 |
| Cultura M 230 V (with a UK plug) | 133192 |
| Cultura M 12 V (with a car plug) | 132697 |
| Thermometer, accessory | 129981 / 05724 |
| Multirack, accessory | 139248 |

IT'S SO EASY TO USE!



Dip or pour

Incubate

Interpret

For more detailed information, please consult the instructions for use.

Orion Diagnostics develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission.


Uricult® is a registered trademark of Orion Diagnostics Oy.



Orion Diagnostics Oy, Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland
Tel. +358 10 4261, Fax +358 10 426 2794, www.oriondiagnostics.com

757-05EN, 10/2017. International version.

Organisms on media



Colony densities *Escherichia coli* CFU/ml

