

QuikRead go  
wrCRP+Hb



# QuikRead go<sup>®</sup> wrCRP+Hb

Test before prescribing

**ORION**  
DIAGNOSTICA

# Test before prescribing

The QuikRead go wrCRP+Hb is a unique wide range point-of-care CRP test, which gives simultaneously reliable and quantitative CRP and haemoglobin (Hb) test results within 2 minutes. The wide range CRP test allows the extension of CRP values from conventional use to testing e.g. neonatal sepsis<sup>1</sup> as well as to predict disease severity in COPD<sup>2</sup> and aid in cardiovascular risk assessment<sup>3</sup>. Room temperature storage for unopened kit and small sample volume ease the everyday work of the healthcare professionals. QuikRead go system is a fully automatic and flexible test system and is connectable to most hospital (HIS) and laboratory information (LIS) systems.

## Because...

### ROOM TEMPERATURE STORAGE

- » Easy to store in the healthcare unit
- » Enables immediate use, no need to warm up the reagents

### WIDE CRP MEASURING RANGE

- » Wide measuring range from whole blood is 0.5–200 mg/l and 0.5–120 mg/l with serum or plasma samples
- » Enables using CRP in wide range of indications

### SMALL SAMPLE VOLUME 10 µL

- » Sample collection is easy and fast
- » Offers improved convenience for the patient

### UNIQUE COMBINATION

- » Hb measuring range 50–220 g/l (5.0–22.0 g/dl, 3.1–13.7 mmol/l)
- » Two results from the same sample give more information for treatment decision

### References

1. Hengst JM: The role of C-reactive protein in the evaluation and management of infants with suspected sepsis. *Adv Neonatal Care* 2003; 3:3–13.
2. Mendy A, Forno E, Niyonsenga T, Gasana J. *Clin Respir J*. 2017. Blood biomarkers as predictors of long-term mortality in COPD. <https://doi.org/10.1111/crj.12752>
3. Kaptoge S, Di Angelantonio E, Lowe G, Pepys MB, Thompson SG, Collins R, Danesh J. C-reactive protein concentration and risk of coronary heart disease, stroke, and mortality: an individual participant meta-analysis. *Emerging Risk Factors Collaboration*. 2010 *Lancet* 375(9709):132-40.



**QuikRead go**<sup>®</sup>

## IT'S SO EASY TO USE!



Sample



Measurement



Result

*For more detailed information, please consult the instructions for use.*

### Product

### Cat. no.

QuikRead go wrCRP+Hb Kit, 50 tests	
» With capillaries and plungers	146255
QuikRead go wrCRP Control 1 ml	153038
QuikRead go wrCRP Control High 1 ml	153039
Ultrasensitive CRP Control 1 ml	68257
QuikRead go Hb 10 µl Control 1 ml	153656
QuikRead go Instrument	133893



Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission.

QuikRead go<sup>®</sup> is a registered trademark of Orion Diagnostica Oy.

Orion Diagnostica Oy  
Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland  
[www.oriondiagnostica.com](http://www.oriondiagnostica.com), [www.quikread.com](http://www.quikread.com), [www.tacklearm.com](http://www.tacklearm.com)