



QuikRead go[®] Strep A

Objective detection of Strep A pharyngitis

Objective detection of Strep A pharyngitis

QuikRead go Strep A is a qualitative, immunochemical test for detection of *Streptococcus pyogenes* (Strep A) in a throat swab sample. The test is easy to use and gives a positive result for throat swabs containing Strep A in concentrations higher than 7×10^4 CFU. The accurate and timely result is available during the patient's consultation to support the decision of prescribing antibiotics.

Because...

FAST

- » Results in less than 4 minutes measurement time
- » Clear negative and positive results are displayed in less than 2 minutes
- » Instrument read results make the visual interpretation unnecessary

FLEXIBLE

- » Automated measurement procedure
- » Recording of the results is easy to fit into the laboratory workflow

TRACEABLE RESULTS

- » User and patient ID can be stored with the results

PROVEN QUIKREAD QUALITY

- » Results automatically stored in the instrument's memory
- » System connectivity features increase reliability of results' handling



QuikRead go

IT'S SO EASY TO USE!



Sample preparation



Measurement



Result

For more detailed information, please consult the instructions for use.

Product

Cat. no.

QuikRead go Strep A, 50 tests	135883
QuikRead go Instrument	133893



Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission. QuikRead go® is a registered trademark of Orion Diagnostica Oy.

Orion Diagnostica Oy
Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland
Tel. +358 10 4261, Fax +358 10 426 2794
www.oriondiagnostica.com, www.quikread.com

885-06GB, 03/2016. International version. Not intended for US market.