

QuikRead go wrCRP



QuikRead go[®] wrCRP

Test before prescribing



Test before prescribing

The QuikRead go wrCRP is a wide range point-of-care CRP test, which gives reliable and quantitative test result within 2 minutes. The wide range CRP test allows the extension of CRP values from conventional use to testing e.g neonatal sepsis¹ as well as to predict disease severity in COPD² and aid in cardiovascular risk assessment³. Room temperature storage for unopened kit and small sample volume ease the everyday work of the healthcare professionals. QuikRead go system is a fully automatic and flexible test system and is connectable to most hospital (HIS) and laboratory information (LIS) systems.

Because...

ROOM TEMPERATURE STORAGE

- » Easy to store in the healthcare unit
- » Enables immediate use, no need to warm up the reagents

WIDE CRP MEASURING RANGE

- » Wide measuring range from whole blood is 0.5–200 mg/l and 0.5–120 mg/l with serum or plasma samples
- » Enables using CRP in wide range of indications

SMALL SAMPLE VOLUME 10 µL

- » Sample collection is easy and fast
- » Offers improved convenience for the patient

References

1. Hengst JM: The role of C-reactive protein in the evaluation and management of infants with suspected sepsis. *Adv Neonatal Care* 2003; 3:3–13.
2. Mendy A, Forno E, Niyonsenga T, Gasana J. *Clin Respir J*. 2017. Blood biomarkers as predictors of long-term mortality in COPD. <https://doi.org/10.1111/crj.12752>
3. Kaptoge S, Di Angelantonio E, Lowe G, Pepys MB, Thompson SG, Collins R, Danesh J. C-reactive protein concentration and risk of coronary heart disease, stroke, and mortality: an individual participant meta-analysis. *Emerging Risk Factors Collaboration*. 2010 *Lancet* 375(9709):132-40.



QuikRead go®

IT'S SO EASY TO USE!



Sample



Measurement



Result

For more detailed information, please consult the instructions for use.

Product

Cat. no.

QuikRead go wrCRP Kit, 50 tests	
» With capillaries and plungers	146521
QuikRead go wrCRP, 500 tests	146522
QuikRead go wrCRP Control 1 ml	153038
QuikRead go wrCRP Control High 1 ml	153039
Ultrasensitive CRP Control 1 ml	68257
QuikRead go Instrument	133893



Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission.

QuikRead go® is a registered trademark of Orion Diagnostica Oy.

Orion Diagnostica Oy
Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland
www.oriondiagnostica.com, www.quikread.com, www.tacklearm.com