



Ordering information

Product	Cat. no.
Easicult Combi 10 tests	67987
Easicult Combi PP 10 tests	05984
Easicult M 10 tests	67686
Easicult TTC 10 tests	67683
Easicult TTC PP 10 tests	05988
Cultura M 230 V	129997/05866
Cultura M 110/115 V	132754/06266
Cultura M 230 V (with a UK plug)	133192
Cultura M 12 V (with a car plug)	132697
Thermometer, accessory	129981/05724
Multitrack, accessory	139248

For more information, please consult Instructions for use.

Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission. Easicult® is a registered trademark of Orion Diagnostica Oy.



Orion Diagnostica Oy
 Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland
 Tel. +358 10 4261, Fax +358 10 426 2794
 www.oriondiagnostica.com

8086-02EN, 01/2018. International version.

It's so easy to use!



Dip or swab



Incubate

Interpret



Easicult® Products

Easicult® Combi Easicult® M Easicult® TTC



Easy and handy microbial controlling for various industries

Easicult is a dipslide culture method that is intended for monitoring microbial contaminations in industrial fluids. The test can easily be used on-site without any microbiological expertise and the slides can be used as convenient transport media if further analysis of the samples is needed.

Reliable

- » High quality products with users all over the world
- » Established method with long experience

Time-saving

- » Easy on-site testing
- » No microbiological skills or special training needed
- » Simple results interpretation – no counting of colonies

Flexible

- » Long shelf-life compared to culture plates
- » Transportable
- » Storage at room temperature
- » Incubation at ambient temperature or in incubator for faster results

Easicult for easy monitoring of microbial contaminations

- » Paper and pulp industry fluids
- » Metal industry fluids
- » Industrial cutting fluids
- » Fuels, lubricating oils, diesel, biodiesel
- » Paints, pigments, varnishes
- » Manufacturers of cars and other motor vehicles
- » Cooling towers, cooling waters
- » Cosmetics industry
- » Fuel tanks
 - » Aircraft industry
 - » Ships and boats



Cultura M incubator is for easy and handy culturing on-site with temperature range of 25–45°C.

Thermometer and multitrack for easy incubation of the tubes are also available.

Easicult® Combi

- » TTC agar for total bacterial counts and Rose Bengal agar for yeasts and moulds
- » For easy monitoring of general hygiene
- » Two product versions: Polypropylene (PP) tube for longer stability and polystyrene (PS) tube for varying temperature conditions



Easicult® M

- » Malt agar for detection of yeasts and moulds
- » For monitoring fungal contamination



Easicult® TTC

- » TTC agar for estimation of total bacterial counts
- » Red coloured colonies make the results interpretation easy
- » Two product versions: Polypropylene (PP) tube for longer stability and polystyrene (PS) tube for varying temperature conditions

