


NEW!




SAMPLE COLLECTOR


for QuikRead go[®] easy CRP

Unique one-step Sample Collector


Unique one-step Sample Collector used with the QuikRead go[®] easy CRP minimizes hands-on steps and enables fast and easy testing at point-of-care.



Take the sample.



Insert Sample Collector in to the cuvette.



CRP result in 2 minutes.

Sample Collector is a sampling device developed specially for QuikRead go easy CRP. The QuikRead go easy CRP kit can be stored at room temperature (up to +25°C) which eases everyday work of healthcare professionals.

The QuikRead go easy CRP kit contains Sample Collectors and is used with the QuikRead go multianalyte point-of-care system for quantitative measurement of CRP.



Available in January 2019.

Test before prescribing



www.oriondiagnostica.com



www.quikread.com



www.tackleamr.com