

New!

QuikRead® iFOB Sampling Set

Compatible combinations


QuikRead® iFOB Sampling Set For use with QuikRead go® iFOB kits IVD 2°C 25°C Tubes for dilution and filtration of faecal samples Orion Diagnostics Oy Kotivuoratie 10, E 8 PO-Box 81, FIN-02101 Espoo, Finland Tel. +358 10 4261, Fax +358 10 426 2794 www.oriondiagnostics.com	REF 151052	LOT
	CONT	∇/50
	QuikRead iFOB Sampling vial	50 x 2 ml
	Patient instruction	50
	Label	50
	Danger	CONT CH ₂ O
		CONT NaN



QuikRead® iFOB Sampling vial 2 ml IVD +2°C +25°C CONT CH ₂ O 0.1 % CONT NaN < 0.1 % Orion Diagnostics Oy, Espoo, Finland	LOT
	1-680815



REF 151051





REF 152390

151052 QuikRead iFOB Sampling set can be used in combination with:


Material no.	Product name
151051	QuikRead go iFOB kit
152390	QuikRead go iFOB Control Quantitative
06027	QuikRead FOB Positive Control


QuikRead® FOB Sampling Set

Compatible combinations

QuikRead® FOB Sampling Set For use with QuikRead® FOB quantitative & QuikRead go® iFOBT kits IVD  25°C Tubes for dilution and filtration of faecal samples  Orion Diagnostica Oy Kotusmäntienkatu 6 B P.O.Box 83, FI-02101 Espoo, Finland Tel. +358 10 6261, Fax +358 10 426 2794 www.oriondiagnostica.com	REF 138822	LOT	
	CONT	▽/60	LOT
	QuikRead FOB Sampling vial	50 x 2 ml	
	Patient instruction	50	
	Label	50	
	CONT NaN ₃		



QuikRead® FOB Sampling vial IVD LOT BUF 2 ml  +25°C CONT NaN ₃ < 0.1 % Orion Diagnostica Oy Espoo, Finland	145234-1
--	----------



138822 QuikRead FOB Sampling set can be used in combination with:

Material no.	Product name
136189	QuikRead go iFOBT Kit
138820	QuikRead FOB quantitative Kit
144675	QuikRead FOB Control Quantitative
06027	QuikRead FOB Positive Control