



QuikRead go[®]

Your support in treatment decisions

Your support in treatment decisions

QuikRead go Instrument	Cat.No.
QuikRead go Instrument (FI, SE, NO, DK, CZ, SK, PL, HU)	135867
QuikRead go Instrument (other languages)	133893
Information about QuikRead go accessories: www.oriondiagnostica.com	

QuikRead go is a fully automatic, reliable and easy to use instrument for measuring the QuikRead go reagent kits.

QuikRead go CRP	Cat.No.
QuikRead go CRP, 50 tests including capillaries and plungers	135171
QuikRead CRP Control	68296
QuikRead go CRP Control High	137071

QuikRead go CRP test is intended for quantitative determination of CRP in whole blood, plasma and serum. The measuring range of the CRP test is 5–200 mg/l.

QuikRead go CRP+Hb	Cat.No.
QuikRead go CRP+Hb, 50 tests including capillaries and plungers	140068
QuikRead CRP Control	68296
QuikRead go CRP Control High	137071
QuikRead go Hb Control	141154

QuikRead go CRP+Hb test is intended for quantitative determination of CRP in whole blood, plasma and serum and quantitative determination of Hb in whole blood. The test gives both CRP and Hb results from one whole blood sample during one analysis. The measuring range of CRP is 5–200 mg/l and the measuring range of Hb is 50–245 g/l (3.1–15.2 mmol/l).

QuikRead go Strep A	Cat.No.
QuikRead go Strep A, 50 tests including controls	135883

QuikRead go Strep A test is intended for detection of *Streptococcus pyogenes* in throat samples when Strep A pharyngitis (tonsillitis) is suspected.

QuikRead go iFOBT	Cat.No.
QuikRead go iFOBT, 50 tests	151051
QuikRead iFOB Sampling set, 50 tests	151052
QuikRead go iFOBT Control Quantitative	152390
QuikRead FOB Positive Control	06027

QuikRead go iFOBT test is intended for detection and quantification of Hb originating from occult blood in human faeces. Concentrations above 15 µg/g of faeces are identified as positive and the instrument displays quantitative values in the range of 15–200 µg/g of faeces. The quantitative values can also be obtained as ng/ml of faeces.

Reference:
Instructions for use.



Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission. QuikRead go® is a registered trademark of Orion Diagnostica Oy.

Orion Diagnostica Oy
Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland
Tel. +358 10 4261, Fax +358 10 426 2794
www.oriondiagnostica.com, www.quikread.com

8014-08EN, 12/2017. International version. Not intended for US market.