



QuikRead go[®] CRP+Hb

Get more information for diagnosis

Get more information for diagnosis

QuikRead go CRP+Hb test is a reliable and easy to use test for quantitative measurement of C-reactive protein (CRP) and haemoglobin (Hb). In addition to clinical signs and symptoms, measurement of CRP can effectively guide antibiotic use. Monitoring CRP levels also provides means for assessing treatment response. In addition, Hb measurement gives information about the patient's condition in both acute and basic healthcare settings. Combination of CRP and Hb results can give valuable information to help differentiating evolving infection from anaemia, as these conditions may have somewhat similar symptoms. Combination test may ease the healthcare practitioner to make treatment decisions especially when the patients cannot specify their symptoms accurately.

Because...

ONE ANALYSIS — TWO RESULTS

- » Reliable CRP and Hb results in single analysis with whole blood samples
- » Analysis of two parameters with one instrument
- » CRP measuring range 5–200 mg/l and Hemoglobin 50–245 g/l

PATIENT FRIENDLY

- » Only one sample is needed, two results available
- » The need of one sample minimizes patient discomfort

EASY TO USE

- » Easy testing procedure
- » Automatic hematocrit correction for CRP results



QuikRead go[®]

IT'S SO EASY TO USE!



Sample



Measurement



Result

For more detailed information, please consult the instructions for use.

Product	Cat. no.
QuikRead go CRP+Hb Kit, 50 tests	
» With capillaries and plungers	140068
QuikRead CRP Control 1 ml	68296
QuikRead go CRP Control High 1 ml	137071
QuikRead go Hb Control	141154
QuikRead Capillaries, 50 pcs	67962
QuikRead Plungers, 50 pcs	67966
QuikRead go Instrument	133893



Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission.

QuikRead go[®] is a registered trademark of Orion Diagnostica Oy.

Orion Diagnostica Oy
Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland
Tel. +358 10 4261, Fax +358 10 426 2794
www.oriondiagnostica.com, www.quikread.com

897-04GB, 03/2016. International version. Not intended for US market.