

Effective Date 03.07.2015  
Revision 2.0

**SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING**

**1.1. Product identifier**

Trade name: UniQ PIIINP RIA  
Catalogue number: 68570

**1.2. Relevant identified uses of the substance or mixture and uses advised against**

The uses of the chemical: For in vitro diagnostic use

**1.3. Details of the supplier of the safety data sheet:**

Manufacturer/importer: Orion Diagnostica Oy  
Street address: Koivu-Mankkaan tie 6 B  
Post-office box: P.O. Box 83  
Postcode: 02101 Espoo, Finland  
Telephone number: +358 10 4261  
Telefax: +35810 426 2794  
E-mail address: product.support@oriondiagnostica.fi  
VAT Reg. No: FI18552161

**1.4. Emergency telephone number**

Orion Diagnostica Oy +358 10 426 3344

**SECTION 2: HAZARDS IDENTIFICATION**

**2.1. Classification of the substance or mixture**

Effective Date 03.07.2015  
 Revision 2.0

Name of the Component	Classification according to regulation EC 1272/2008	H Statements
PIIINP RIA tracer	Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008. ;	
PIIINP RIA antiserum	Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008. ;	
PIIINP RIA calibrators	Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008. ;	
PIIINP RIA control	Aquatic Chronic 3 ;H412	Harmful to aquatic life with long lasting effects.
Procollagen separation reagent	Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008. ;	

## 2.2. Label elements

Name of the Component	Hazard Pictograms	Signal word	H Statements	P Statements
PIIINP RIA tracer		Not applicable	NA Not applicable	NA Not applicable
PIIINP RIA antiserum		Not applicable	NA Not applicable	NA Not applicable
PIIINP RIA calibrators		Not applicable	NA Not applicable	NA Not applicable
PIIINP RIA control		Not applicable	H412 Harmful to aquatic life with long lasting	P273 Avoid release to the environment.

Effective Date 03.07.2015  
 Revision 2.0

			effects.	P501 Dispose of contents according to national and local law.
Procollagen separation reagent		Not applicable	NA Not applicable	NA Not applicable

**2.3. Other hazards**

-

**SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS**

**3.2. Mixtures**

Name of the ingredient	CAS-number:	EC-number	Reach registration number:	Concentration/ Limit	Classification according to regulation EC 1272/2008
------------------------	-------------	-----------	----------------------------	----------------------	---

**PIIINP RIA antiserum**

Sodium azide	26628-22-8		NA	<0,1%	Acute Tox. 2, Aquatic Acute 1, Aquatic Chronic 1 ;H300 EUH032 H400 H410
--------------	------------	--	----	-------	---

**PIIINP RIA calibrators**

Sodium azide	26628-22-8		NA	<0,1%	Acute Tox. 2, Aquatic Acute 1, Aquatic Chronic 1 ;H300 EUH032 H400 H410
--------------	------------	--	----	-------	---

Effective Date 03.07.2015  
 Revision 2.0

**PIIINP RIA control**

Sodium azide	26628-22-8		NA	<1%	Acute Tox. 2, Aquatic Acute 1, Aquatic Chronic 1 ;H300 EUH032 H400 H410
--------------	------------	--	----	-----	--

**PIIINP RIA tracer**

Sodium azide	26628-22-8		NA	<0,1%	Acute Tox. 2, Aquatic Acute 1, Aquatic Chronic 1 ;H300 EUH032 H400 H410
--------------	------------	--	----	-------	--

**Procollagen separation reagent**

Sodium azide	26628-22-8		NA	<0,1%	Acute Tox. 2, Aquatic Acute 1, Aquatic Chronic 1 ;H300 EUH032 H400 H410
--------------	------------	--	----	-------	--

**Other information**

Rehydrated reagent contains < 0,1% sodium azide.

**SECTION 4: FIRST AID MEASURES**

**4.1. Description of first aid measures**

**Skin contact**

Wash with soap and water.

**Eye contact**

Rinse continuously with water for several minutes.

Effective Date 03.07.2015  
Revision 2.0

**Ingestion** If the patient is conscious, give water (up to 2 glasses).  
Give charcoal.

#### 4.2. Most important symptoms and effects, both acute and delayed

Data not available.

#### 4.3. Indication of any immediate medical attention and special treatment needed

Data not available.

### SECTION 5: FIRE-FIGHTING MEASURES

#### 5.1. Extinguishing media

**Suitable extinguishing media** Choose suitable extinguishing media according to the environment.

#### 5.2. Special hazards arising from the substance or mixture

Possibility to small amounts of harmful gases or vapours.

#### 5.3. Advice for fire-fighters

**Special protective equipment for fire-fighters** No special protective equipment needed.

### SECTION 6: ACCIDENTAL RELEASE MEASURES

#### 6.1. Personal precautions, protective equipment and emergency procedures

Wear protective clothing.

#### 6.2. Environmental precautions

No special environmental precautions needed.

#### 6.3. Methods and material for containment and cleaning up

In case of spillage clean with paper towel and disinfect.

Effective Date 03.07.2015  
Revision 2.0

#### 6.4. Reference to other sections

See Section 8 and 13

### SECTION 7: HANDLING AND STORAGE

#### 7.1. Precautions for safe handling

Do not eat, drink or smoke at workplace. Wash hands after working with substance.

#### 7.2. Conditions for safe storage, including any incompatibilities

Store at +2 - 8 °C.

#### 7.3. Specific end use(s)

No information identified.

### SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

#### 8.1. Control parameters

International OEL values Sodium azide 0,1 mg/m<sup>3</sup> (8h) , 0,3 mg/m<sup>3</sup> (15 min)

#### 8.2. Exposure controls

##### Hand protection

Wear protective gloves.

##### Skin protection

Wear protective clothing.

### SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

#### Appearance

Physical state	Liquid / Solid
Colour	Antiserum blue Tracer red Other reagents colourless

Effective Date 03.07.2015  
Revision 2.0

**Odour** Odourless

**9.1. Information on basic physical and chemical properties**

**Water solubility** Reagents soluble

**9.2. Other information**

-

**SECTION 10: STABILITY AND REACTIVITY**

**10.1. Reactivity**

-

**10.2. Chemical stability**

-

**10.3. Possibility of hazardous reactions**

-

**10.4. Conditions to avoid**

-

**10.5. Incompatible materials**

When in contact with acids, sodium azide forms highly toxic gas. Sodium azide reacts with metals to form explosive metal azides.

Effective Date 03.07.2015  
Revision 2.0

## 10.6 Hazardous decomposition products

-

## SECTION 11: TOXICOLOGICAL INFORMATION

### 11.1. Information on toxicological effects

#### Acute toxicity

Quantitative data on the toxicity for this product are not available.

Sodium azide: LD50 oral (rat) 27 mg/kg, LD50 dermal (rabbit) 20 mg/kg

#### Skin irritation and corrosion

Data not available.

#### Skin sensitization

Data not available.

#### Serious eye damage and irritation

Data not available.

#### Respiratory irritation

Data not available.

#### Respiratory sensitization

Data not available.

#### Carcinogenicity

Data not available.

#### Germ cell mutagenicity

Data not available.

#### Reproductive toxicity

Data not available.

#### STOT-single exposure

Data not available.

#### STOT-repeated exposure

Data not available.

#### Aspiration hazard

Data not available.

#### Repeated dose toxicity

Data not available.



Effective Date 03.07.2015  
Revision 2.0

## SECTION 12: ECOLOGICAL INFORMATION

### 12.1. Toxicity

#### Aquatic toxicity

Quantitative data on the toxicity for this product are not available.  
Fish toxicity of sodium azide: L. Macrochirus LC50 0,7 mg/l/96 h. Daphnia toxicity of sodium azide : Daphnia pulex EC50 4,2 mg/l/48 h.

#### Toxic effects on other organisms

Data not available.

### 12.2. Persistence and degradability

Data not available.

### 12.3. Bioaccumulative potential

Data not available.

### 12.4. Mobility in soil

Data not available.

### 12.5. Results of PBT and vPvB assessment

Data not available.

### 12.6. Other adverse effects

Data not available.

## SECTION 13: DISPOSAL CONSIDERATIONS

### 13.1. Waste treatment methods

Used product should be considered infectious and should be handled respectively.  
Disposal of all sample and test material should be done in compliance with national, state and local regulations.  
If not officially differently specified, packaging may be treated like household waste or recycled.

Effective Date 03.07.2015  
Revision 2.0

**SECTION 14: TRANSPORT INFORMATION****14.1. UN Number**

UN 2910

**14.2. UN Proper shipping name**

Radioactive material, excepted package - limited quantity of material ,

**14.3. Transport hazard class(es) ( ADR/RID,IMDG,ICAO/IATA)**

-

**14.4. Packing group**

-

**14.5. Environmental hazard**

-

**14.6. Special precautions for user**

-

**14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC code**

Not applicable

**SECTION 15: REGULATORY INFORMATION****15.1. Safety, health and environmental regulations/legislation specific for substance or mixture**

-

**15.2. Chemical Safety Assessment**

No

**SECTION 16: OTHER INFORMATION****List of H statements**

H300, EUH032, H400, H410      Fatal if swallowed.  
Contact with acids liberates very toxic gas.  
Very toxic to aquatic life.

Effective Date 03.07.2015  
Revision 2.0

Very toxic to aquatic life with long lasting effects.

H300, EUH032, H400, H410 Fatal if swallowed.  
Contact with acids liberates very toxic gas.  
Very toxic to aquatic life.  
Very toxic to aquatic life with long lasting effects.

H300, EUH032, H400, H410 Fatal if swallowed.  
Contact with acids liberates very toxic gas.  
Very toxic to aquatic life.  
Very toxic to aquatic life with long lasting effects.

H300, EUH032, H400, H410 Fatal if swallowed.  
Contact with acids liberates very toxic gas.  
Very toxic to aquatic life.  
Very toxic to aquatic life with long lasting effects.

H300, EUH032, H400, H410 Fatal if swallowed.  
Contact with acids liberates very toxic gas.  
Very toxic to aquatic life.  
Very toxic to aquatic life with long lasting effects.

#### Training advice

Read Instructions for Use

#### Further information

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. Orion Diagnostica Oy shall not be held liable for any damage resulting from handling or from contact with the above product.

#### Sources of key data used to compile the Safety Data Sheet

Directive 1272/2008/EC.  
SDS for Sodium azide.  
Instructions for use.

#### Information which has been added, deleted or revised

Updated to meet the CLP Regulation