



QuikRead go[®] CRP

Test before prescribing

Test before prescribing

QuikRead go CRP is an easy-to-use test for quantitative measurement of C-reactive protein (CRP). CRP is normally present at very low concentrations in the blood of healthy people¹. CRP concentrations are markedly increased in bacterial infections, whereas viral infections mostly induce a very modest elevation or none at all². In addition to clinical signs and symptoms, measurement of CRP can therefore effectively guide antibiotic use³. Monitoring CRP levels also provides an objective means of assessing treatment response, as CRP levels fall rapidly as a result of effective therapy².

Because...

RELIABLE

- » Proven quality – comparable to laboratory test results⁴
- » Reproducible quantitative CRP result within range of 5–200 mg/l⁵

QUICK

- » Fast point-of-care CRP system⁴
- » Test time 2 minutes

EASY TO USE

- » Easy testing procedure⁴
- » Automatic hematocrit correction
- » Performed on a finger prick blood sample (alternatively venous blood, plasma or serum)

References

1. Shine B et al. Solid phase radioimmunoassays for human C-reactive protein. Clin Chim Acta. 1981;117(1):13-23.
2. Pepys, M.B., The Acute Phase response and C-reactive Protein, in Oxford Textbook of Medicine, D.A. Warrell, et al., Editors. 2003, Oxford University Press. p. 150-156.
3. Aabenhus, R., et al., Biomarkers as point-of-care tests to guide prescription of antibiotics in patients with acute respiratory infections in primary care. Cochrane Database Syst Rev. 2014;11.
4. Data on file, Orion Diagnostica.
5. QuikRead go CRP Instructions for use.



QuikRead go[®]

IT'S SO EASY TO USE!



Sample



Measurement



Result

For more detailed information, please consult the instructions for use.

Product

Cat. no.

QuikRead go CRP Kit, 50 tests	
» With capillaries and plungers	135171
QuikRead CRP Control 1 ml	68296
QuikRead go CRP Control High 1 ml	137071
QuikRead Capillaries, 50 pcs	67962
QuikRead Plungers, 50 pcs	67966
QuikRead go Instrument	133893



Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission.

QuikRead go[®] is a registered trademark of Orion Diagnostica Oy.

Orion Diagnostica Oy
Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland
Tel. +358 10 4261, Fax +358 10 426 2794
www.oriondiagnostica.com, www.quikread.com

862-05GB, 03/2017. International version. Not intended for US market.